



Orthopaedics

High power efficiently clears bone and cement debris during total joint procedures. Specific tips are available for hip and knee applications.

Trauma / Wound Care

Variety of tips provide a wide range of pressure options for soft tissue treatment. Low power setting provides gentle lavage to help remove necrotic tissue and facilitate healing.



The *Pulsavac Plus* System offers a disposable alternative to the proven *Pulsavac*® Wound Debridement System.

KitPacks

- 5150-420 Component Kit
- 5150-475 Fan Spray Kit
- 5150-482 Hip Kit
- 5150-495 ShowerSprayKit

Multi-Orifice Tips

- † 5150-175 High Capacity Fan Spray Tip w/Splash Shield
- †† 5150-176 High Capacity Fan Spray Tip w/Filter
- † 5150-195 High Capacity Shower Spray Kit w/Splash Shield, 12.7cm length

Debridement / Suction Tips

- † 5150-155 High Capacity Single Stream Tip w/Splash Shield, 12.7cm length
- 5150-065 Large Fluid Orifice Tip w/Wide Single Stream, 7.2cm length

Femoral Tips

- 5150-181 High Capacity IM Tip Four Orifices, 20.32cm x 1.27cm
- 5150-182 High Capacity IM Tip w/Radial Spray, 22.86cm x .89cm
- 5150-183 Capacity IM Brush Tip w/Radial Spray, 22.86cm x .89cm x 1.52cm
- †† 5150-184 High Capacity IM Hip Revision Tip, 23cm length
- 5150-185 High Capacity IM Tip Four Orifices, 30.48cm x 1.27cm

Specialty Tip

- 5150-145 High Capacity Adjustable Length Tip w/Splash Shield

† U.S. Patent 4,692,140
 †† U.S. Patents D365,146; 5,499,970

For more information, contact your Zimmer representative.

CE UL This product is latex-free.



www.zimmer.com



97-5160-11 Rev. 1 20VL - Printed in U.S.A. ©2000 Zimmer, Inc.