

# WALKERS

## ZIMMER® FIXED HINGE LOWER LEG WALKER

- Contoured strut design allows the walker frame to conform to patient's anatomy ensuring a better fit
- Shock absorbing sole reduces impact from the heel strike to improve comfort during ambulation
- Rocker bottom for plantar unloading is engineered to be low profile and wide to minimize leg length discrepancy and help promote a natural gait with increased stability
- Plastic uprights provide support and increased protection
- Walker provides support and comfort during ambulation and replaces a walking cast
- Indicated for soft tissue injuries (grade II and III sprains), stable fractures, post-operative use, trauma, rehabilitation and other conditions recommended by a physician
- Lightweight design is easier on patients



### WALKER

Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-004-00	Small	6 - 8	4 1/2 - 7
00-8916-005-00	Medium	8 1/2 - 11 1/2	7 1/2 - 10 1/2
00-8916-006-00	Large	11 1/2 - 13 1/2	10 1/2 - 12 1/2

### REPLACEMENT LINERS

Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-004-01	Small	6 - 8	4 1/2 - 7
00-8916-005-01	Medium	8 1/2 - 11 1/2	7 1/2 - 10 1/2
00-8916-006-01	Large	11 1/2 - 13 1/2	10 1/2 - 12 1/2

## ZIMMER® LOW TOP ANKLE WALKER

- Lowered and contoured strut for when a taller version is not needed
- Rocker bottom for plantar unloading is engineered to be low profile and wide to minimize leg length discrepancy and help promote a natural gait with increased stability



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-041-00	Small	6 - 8	4 1/2 - 7
00-8916-042-00	Medium	8 1/2 - 11 1/2	7 1/2 - 10 1/2
00-8916-043-00	Large	11 1/2 - 13 1/2	10 1/2 - 12 1/2



# WALKERS

## ZIMMER® AIR WALKER

- Pneumatic system increases stability while decreasing pain and swelling
- Easy to use pump for inflation/deflation
- Air is easily adjusted to achieve desired compression and custom fit for better patient compliance
- Accommodates the swelling patterns that occur throughout the rehabilitation process as well as dressings
- Rocker bottom is specially engineered to be low and wide to help promote a natural gait and increased stability



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-024-00	Small	6 - 8	4 1/2 - 7
00-8916-025-00	Medium	8 1/2 - 11 1/2	7 1/2 - 10 1/2
00-8916-026-00	Large	11 1/2 - 13 1/2	10 1/2 - 12 1/2

## ZIMMER® LOWER LEG ORTHOSIS II

- Hinged ankle design provides up to 45 degrees of controlled ankle motion
- Low profile design provides better anatomical conformity for smoother gait
- Replacement liners are available



WALKER			
Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-014-00	Small	6 1/2 - 9 1/2	Up to 7
00-8916-015-00	Medium	10 - 12 1/2	8 - 9
00-8916-016-00	Large	13 - 15	10 and up

REPLACEMENT LINERS			
Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-011-00	Small	6 1/2 - 9 1/2	Up to 7
00-8916-012-00	Medium	10 - 12 1/2	8 - 9
00-8916-013-00	Large	13 - 15	10 and up

## DH PRESSURE RELIEF™\* WALKER

- Pre-inflated bladder for strut protection
- 90 degrees fixation of ankle to reduce shearing
- Selective pressure relief easily created by removal of independent shock-absorbing hexagons
- Rocker bottom is specially engineered to be low and wide to help promote a natural gait and increased stability
- Effective solution for treating mild to moderate plantar ulcers
- Inner bootie designed to consistently locate the foot relative to the pressure relief insole
- Patient-friendly leg straps
- Liner and sole are easily removed for cleaning



Prod. No.	Size	Shoe Size	
		Women's	Men's
00-8916-031-00	Small	Up to 10	Up to 8
00-8916-032-00	Medium	10 1/2 - 12 1/2	8 1/2 - 10 1/2
00-8916-033-00	Large	12 1/2 and up	10 1/2 and up

\*Trademark of Royce Medical Company